



HOLON APP QUICK START GUIDE


Version 02
July 2024

→] Launching the Holon App & Signing In




Click the Holon icon on your desktop to open the app.

2



On the Profiles page, select **New User**.


3



Enter your **User ID** & **Password** and click **Sign in**.

🌐 Sensing Against Your EHR


4



Open a patient chart in your EHR to get started. Holon will automatically sense the patient and load data into the app.

🔔 Notifications

5




Once data is ready in the app, a notification will appear on your screen detailing the patient's available data.

- Click a category (e.g. **Care Gaps**, **Coding Gaps**, **Patient Summary**) to display content
- Click **Close** to dismiss
- Red dots on icons indicate data is available


📍 Navigation

6




Click category icons under **Patient Care** to access content.

7



Use these icons to control the app to your preferences.

8



Keyboard shortcut "ctrl+shift+M" will reset the app for dual monitors or if it gets dragged off screen.

🔗 Content

Category content comes from your source data system. Category icons and names are configured by your admin and may differ from this Quick Start Guide. If you have any questions on data content, please contact your admin for in depth training material.

❓ Support

If you have any questions or need app support, visit

help.holonsolutions.com

Log in with your Holon app credentials.